



"A new hope for a better tomorrow"

HILZO STENTS CATALOGUE



HILZO
s t e n t s



Colonic Stents



Colonic Stents

Used as a palliative treatment for patients with colonic strictures.



1 Uncovered type

- Special cell construction provides excellent flexibility for placements in colonic areas.

Art. No.	Stent		Introducer System	
	Diameter (mm)	Length (cm)	Diameter	Usable Length (cm)
CV2208	22	8	10Fr	220
CV2210		10		
CV2212		12		
CV2408	24	8		
CV2410		10		
CV2412		12		

! 20, 26, 28mm body diameter available upon doctor's request.

Colonic Stents

Used as a palliative treatment for patients with colonic strictures.



2 Covered type

- PTFE covered type is designed to prevent tumor ingrowth (longer patency).
- Longer diameter heads are designed to reduce migration.

Fully covered (PTFE body + Silicone covered ends)

Art. No.	Stent					Introducer System	
	Diameter (mm)	Length (cm)	Head			Diameter	Usable Length (cm)
			Diameter (mm)	Length (cm)			
				L	R		
CVP2008HC1	20	8	28	1.5	1.5	10.5Fr	220
CVP2010HC1		10					
CVP2012HC1		12					

! · 18mm body diameter available upon doctor's request.

Partially covered (PTFE body + uncovered ends)

Art. No.	Stent					Introducer System	
	Diameter (mm)	Length (cm)	Head			Diameter	Usable Length (cm)
			Diameter (mm)	Length (cm)			
				L	R		
CVP2008HU1	20	8	28	1.5	1.5	10.5Fr	220
CVP2010HU1		10					
CVP2012HU1		12					

! · 18mm body diameter available upon doctor's request.



CONTACT US

sales@bcmglobal.co.kr

TEL. +82-31-926-2727

