



"A new hope for a better tomorrow"

HILZO STENTS CATALOGUE



HILZO
s t e n t s



A close-up photograph of a gloved hand holding a blue and white mesh stent. The stent is a cylindrical, woven mesh structure with a blue outer layer and a white inner layer. The hand is wearing a white glove. The background is dark, and the stent is reflected on a surface below it.

Pyloric / Duodenal Stents



Pyloric/Duodenal Stents

Used as a palliative treatment for patients with pyloric and duodenal strictures.



1 Uncovered type

- Unique cell construction provides good flexibility for placements in tortuous anatomies.
- Longer diameter heads are designed to reduce migration.

| Art. No. | Stent | | | | | Introducer System | |
|----------|---------------|-------------|---------------|-------------|-----|-------------------|--------------------|
| | Diameter (mm) | Length (cm) | Head | | | Diameter | Usable Length (cm) |
| | | | Diameter (mm) | Length (cm) | | | |
| | | | | L | R | | |
| PV2008H | 20 | 8 | 28 | 1.5 | 1.5 | 10Fr | 180 |
| PV2010H | | 10 | | | | | |
| PV2012H | | 12 | | | | | |

! 18mm body diameter available upon doctor's request.

Pyloric/Duodenal Stents

Used as a palliative treatment for patients with pyloric and duodenal strictures.



2 Covered type

- PTFE covered type is designed to prevent tumor ingrowth.
- Longer diameter heads are designed to reduce migration.

Fully covered (PTFE body + Silicone covered ends)

| Art. No. | Stent | | | | | Introducer System | |
|------------|---------------|-------------|---------------|-------------|-----|-------------------|--------------------|
| | Diameter (mm) | Length (cm) | Head | | | Diameter | Usable Length (cm) |
| | | | Diameter (mm) | Length (cm) | | | |
| L | R | | | | | | |
| PVP2008HC1 | 20 | 8 | 28 | 1.5 | 1.5 | 10.5Fr | 180 |
| PVP2010HC1 | | 10 | | | | | |
| PVP2012HC1 | | 12 | | | | | |

Partially covered (PTFE body + uncovered ends)

| Art. No. | Stent | | | | | Introducer System | |
|------------|---------------|-------------|---------------|-------------|-----|-------------------|--------------------|
| | Diameter (mm) | Length (cm) | Head | | | Diameter | Usable Length (cm) |
| | | | Diameter (mm) | Length (cm) | | | |
| L | R | | | | | | |
| PVP2008HU1 | 20 | 8 | 28 | 1.5 | 1.5 | 10.5Fr | 180 |
| PVP2010HU1 | | 10 | | | | | |
| PVP2012HU1 | | 12 | | | | | |

! 18mm body diameter available upon doctor's request.



CONTACT US

sales@bcmglobal.co.kr

TEL. +82-31-926-2727

