

**medax**

**MED-I**

**BONE MARROW ASPIRATION NEEDLE  
WITH ERGONOMIC HANDLE**



**EXPLANT/TRANSPLANT VERSION**

Bone Marrow aspiration needle with ergonomic handle “XI”, also available in the explant/transplant version “EI”.

- Ergonomic handle for a comfortable procedure and a firm grip.
- Triple sharpened tip cannula to allow a non traumatic, fast and easy penetration into the bone marrow cavity.
- Adjustable depth setter.
- Luer Lock connection for syringe.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

**Your partner in Biopsy and Special Needles**

## BONE MARROW ASPIRATION NEEDLE WITH ERGONOMIC HANDLE

### Ordering information:

GAUGE	COLOUR	LENGTH (MM)					
		14 030	14 040	14 050	14 070	14 090	14 110
14	● Green	14 030	14 040	14 050	14 070	14 090	14 110
15	● Blue	15 030	15 040	15 050	15 070	15 090	15 110
16	○ White	16 030	16 040	16 050	16 070	16 090	16 110
18	● Pink	18 030	18 040	18 050	18 070	18 090	18 110

Different sizes on demand

### Ordering information:

MED-I ASPIRATION TIP

X	I	GAUGE		LENGTH		0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

MED-I EXPLANT/TRANSPLANT TIP

E	I	GAUGE		LENGTH		0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

Single box: 10 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
 Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485:2016 standards.  
 In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.  
 All products undergo intensive clinical testing and are fully EC and FDA approved.