

medax

UNICUT

DISPOSABLE BIOPSY NEEDLE
FOR AUTOMATIC REUSABLE BIOPSY SYSTEMS



Disposable biopsy needle compatible with:

BIP EvoCore®

Uromed Corazor®

Bard® Magnum®

- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Trocar tip.
- Can be used with CU introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

UNICUT

DISPOSABLE BIOPSY NEEDLE
FOR AUTOMATIC REUSABLE BIOPSY SYSTEMS

GAUGE	COLOUR	LENGTH (MM)							
12	Light blu	12 080	12 100	12 120	12 160	12 180	12 200	12 250	12 300
14	Green	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300
16	White	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300
18	Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300
20	Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300

Different sizes on demand

Ordering information:

UNCUT NORMAL TIP	U	T	GAUGE	LENGTH	0	0
UNICUT TROCAR TIP	U	T	GAUGE	LENGTH	0	1
INTRODUCER COAXIAL NEEDLE FOR UNICUT	C	U	GAUGE	LENGTH	0	0
UNICUT "NORMAL TIP" "WITH COAXIAL NEEDLE IN THE SAME PACKAGE	U	T	GAUGE	LENGTH	C	0

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

